

**Supplementary Table S1.** Primers used for PCR detection of *Toxoplasma gondii* oocysts in cat fecal samples

Target gene	Primer sequence	Amplicon size
repetitive 529 bp	5'-CGCTGCAGGGAGGAAGACGAAAGTTG-3'	529 bp
	5'-CGCTGCAGACACAGTGCATCTGGATT-3'	
B1 gene	5'-GGAAGTGCATCCGTTTCATGAG-3'	194 bp
	5'-TCTTTAAAGCGTTCGTGGTC-3'	